

Recombinant Human CRIF1 protein ab113863

画像数 1

製品の詳細

製品名	Recombinant Human CRIF1 protein
精製度	> 90 % SDS-PAGE. ab113863 was purified by using conventional chromatography techniques.
発現系	Escherichia coli
アクセッション番号	<u>Q8TAE8</u>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHSSGLVPRGSHMPRWQLGPRYAAKQFARYGA ASGVVPGSLW PSPEQLRELEAEEREWYPSLATMQESLRVKQLAEEQRRERE QHIAECMA KMPQMIVNWQQQRENWEKAQADKERRARLQAEAQELLGYQV DPRSARFQ ELLQDLEKKERKRLKEEKQKRKKEARAAALAAAVAQDPAASG APSS
予測される分子量	23 kDa including tags
領域	48 to 222
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab113863** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Mass Spectrometry
質量分析	MALDI-TOF
製品の状態	Liquid
備考	Previously labelled as GADD45GIP1.

前処理および保存

保存方法および安定性

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.03% DTT, 0.32% Tris HCl, 40% Glycerol (glycerin, glycerine), 1.17% Sodium chloride

関連情報

機能

Acts as a negative regulator of G1 to S cell cycle phase progression by inhibiting cyclin-dependent kinases. Inhibitory effects are additive with GADD45 proteins but occurs also in the absence of GADD45 proteins. Acts as a repressor of the orphan nuclear receptor NR4A1 by inhibiting AB domain-mediated transcriptional activity. May be involved in the hormone-mediated regulation of NR4A1 transcriptional activity.


組織特異性

Widely expressed. Highly expressed in the thyroid gland, heart, lymph nodes, trachea and adrenal tissues. Expressed at lower level in liver skeletal muscle, kidney, pancreas, testis, ovary and stomach. Barely detectable in adrenal adenoma and papillary thyroid cancer.

細胞内局在

Nucleus.

画像



ab113863 at 3 µg analysed by 15% SDS-PAGE.

SDS-PAGE - Recombinant Human CRIF1 protein
(ab113863)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors