

Product datasheet

Recombinant Human CXCR3 protein ab114323

画像数 1

製品の詳細

製品名	Recombinant Human CXCR3 protein
発現系	Wheat germ
アクセッション番号	P49682
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	SGLCKVAGALFNINFYAGALLLACISFDRYLNIVHATQLYRR GPPARVTL TCLAVWGLCLLFALPDFIFLSAHHDERLNATHCQYNFPQVGR TALRVLQL
予測される分子量	37 kDa
領域	121 to 220
タグ	GST tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab114323** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Western blot ELISA
製品の状態	Liquid


前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
------------	---

関連情報

機能	Receptor for CXCL9, CXCL10 and CXCL11 and mediates the proliferation of human mesangial cells (HMC). Isoform 2 is a receptor for CXCL4 and also mediates the inhibitory activities of CXCL9, CXCL10 and CXCL11 on the growth of human microvascular endothelial cells (HMVEC). May play a role in angiogenesis. Isoform 3 mediates the activity of CXCL11.
組織特異性	Isoform 1 and isoform 2 are mainly expressed in heart, kidney, liver and skeletal muscle. Isoform 1 is also expressed in placenta.
配列類似性	Belongs to the G-protein coupled receptor 1 family.
細胞内局在	Cell membrane.

画像



SDS-PAGE analysis of ab114323 on a 12.5% gel stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors