

Recombinant Human DR1 protein ab98246

画像数 1

製品の詳細

| | |
|-------------|---|
| 製品名 | Recombinant Human DR1 protein |
| 精製度 | >= 85 % SDS-PAGE. ab98246 was purified using conventional chromatography. |
| 発現系 | Escherichia coli |
| アクセッション番号 | <u>Q01658</u> |
| タンパク質長 | Full length protein |
| Animal free | No |
| 由来 | Recombinant |
| 生物種 | Human |
| 配列 | <p>MGSSHHHHHHSSGLVPRGSHMASSSGNDDDLTIPRAAINK</p> <p>MIKETLPNVR</p> <p>VANDARELVVNCCTEFIHLISSEANEICNKSEKKTISPEHVI</p> <p>QALESLGF</p> <p>GSYISEVKEVLQECKTVALKRRKASSRLENLGIPEEELLRQQ</p> <p>QELFAKAR</p> <p>QQQAELAQQEWLQMQQAAQQAQLAAASASASNQAGSSQDEED</p> <p>DDDI</p> |
| 予測される分子量 | 22 kDa including tags |
| 領域 | 1 to 176 |
| タグ | His tag N-Terminus |

特性

Our **Abpromise guarantee** covers the use of **ab98246** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|----------|-------------------|
| アプリケーション | SDS-PAGE |
| | Mass Spectrometry |
| 質量分析 | MALDI-TOF |
| 製品の状態 | Liquid |

前処理および保存

保存方法および安定性

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 8.00

Constituents: 0.00174% PMSF, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

関連情報

機能

The association of the DR1/DRAP1 heterodimer with TBP results in a functional repression of both activated and basal transcription of class II genes. This interaction precludes the formation of a transcription-competent complex by inhibiting the association of TFIIA and/or TFIIB with TBP. Can bind to DNA on its own. Component of the ATAC complex, a complex with histone acetyltransferase activity on histones H3 and H4.

配列類似性

Belongs to the NC2 beta/DR1 family.
Contains 1 histone-fold domain.


翻訳後修飾

Phosphorylation regulates its interaction with TBP. Not phosphorylated when bound to DRAP1.

細胞内局在

Nucleus.

画像



3ug by SDS-PAGE under reducing conditions and visualized by coomassie blue stain.

SDS-PAGE - Recombinant Human DR1 protein
(ab98246)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors