

### Recombinant Human E2F3 protein ab114534

画像数 1

#### 製品の詳細

製品名	Recombinant Human E2F3 protein
発現系	Wheat germ
アクセッション番号	<b><u>O00716</u></b>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MQSGGGVKTDSTLNSLCGYAWVYVWEEKQRCRLSSFFSSS ASIPGLLP SHTLDLVQNVGVVLDEALGWGRERELCVKCLLEMHCGVFSCM GNHLCQAF PHFPYLSHLVSCLCFQLCVILFASCTKLIFSKV
予測される分子量	41 kDa including tags
領域	1 to 133

#### 特性

Our **Abpromise guarantee** covers the use of **ab114534** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA
	SDS-PAGE
	Western blot
製品の状態	Liquid


#### 前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
------------	---

#### 関連情報

<b>機能</b>	Transcription activator that binds DNA cooperatively with DP proteins through the E2 recognition site, 5'-TTTC[CG]CGC-3' found in the promoter region of a number of genes whose products are involved in cell cycle regulation or in DNA replication. The DRTF1/E2F complex functions in the control of cell-cycle progression from G1 to S phase. E2F-3 binds specifically to RB1 protein, in a cell-cycle dependent manner.
<b>配列類似性</b>	Belongs to the E2F/DP family.
<b>細胞内局在</b>	Nucleus.

## 画像



12.5% SDS-PAGE showing ab114534 at approximately 40.74kDa stained with Coomassie Blue.

SDS-PAGE - Recombinant Human E2F3 protein (ab114534)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors