

### Recombinant Human Eotaxin 2 protein ab9834

#### 製品の詳細

製品名	Recombinant Human Eotaxin 2 protein
生理活性	Determined by its ability to chemoattract using human peripheral blood eosinophils using a concentration of 50.0-100.0 ng/ml.
精製度	> 98 % SDS-PAGE. Sterile filtered Greater than 98% pure by HPLC analyses. Endotoxin level is less than 0.1 ng per g (1EU/g).
エンドトキシン・レベル	< 0.100 Eu/μg
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	VVIPSPCCMF FVSKRIPENR VVSYQLSSRS TCLKGGVIFT TKKGQFCGD PKQEWVQRYM KNLDAKQKKA SPRARAVA
予測される分子量	13 kDa

#### 特性

Our **Abpromise guarantee** covers the use of **ab9834** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

製品の状態	Lyophilized
備考	The biological activity of this product is determined by its ability to chemoattract human peripheral blood eosinophils using a concentration range of 50.0-100.0 ng/ml.

#### 前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
再構成	For lot specific reconstitution information please contact our Scientific Support Team.

## 関連情報

---

機能	Chemotactic for resting T-lymphocytes, and eosinophils. Has lower chemotactic activity for neutrophils but none for monocytes and activated lymphocytes. Is a strong suppressor of colony formation by a multipotential hematopoietic progenitor cell line. Binds to CCR3.
組織特異性	Activated monocytes and activated T lymphocytes.
配列類似性	Belongs to the intercrine beta (chemokine CC) family.
翻訳後修飾	N-glycosylated.
細胞内局在	Secreted.

---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors