

Recombinant Human FABP6 protein ab87381

画像数 1

製品の詳細

製品名	Recombinant Human FABP6 protein
生理活性	ab87381 corresponds to isoform 2, 128 aas.
精製度	> 95 % SDS-PAGE. Purified by using conventional chromatography techniques.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MAFTGKFEME SEKNYDEFMK LLGISSDVIE KAHNFKIVTE VQQDGGQDFTW SQHYGGHTM TNKFTVGKES NIQTMGGKTF KATVQMEGGK LVVNFPNYHQ TSEIVGDKLV EVSTIGGVTY ERVSKRLA

特性

Our **Abpromise guarantee** covers the use of **ab87381** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid
備考	ab87381 corresponds to isoform 2, 128 aas.


前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C. pH: 7.40 Constituents: PBS, 10% Glycerol (glycerin, glycerine)
------------	---

関連情報

機能	Ileal protein which stimulates gastric acid and pepsinogen secretion. Seems to be able to bind to bile salts and bilirubins. Isoform 2 is essential for the survival of colon cancer cells to bile acid-induced apoptosis.
組織特異性	Isoform 2 is expressed in colorectal adenocarcinomas and their adjacent normal mucosa (at protein level). Isoform 1 is expressed in the jejunum, ileum, cecum and ascending colon intestine. Isoform 2 is expressed in the gallbladder, duodenum, jejunum, ileum, cecum, ascending, transverse and descending colon, sigmoid colon and rectum.
配列類似性	Belongs to the calycin superfamily. Fatty-acid binding protein (FABP) family.
ドメイン	Forms a beta-barrel structure that accommodates the hydrophobic ligand in its interior.
細胞内局在	Cytoplasm and Cytoplasm. Localized close to nucleus on the apical side of both normal and neoplastic cells.

画像



ab87381 on 15% SDS-PAGE (3μg)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors