

Recombinant Human FAM19A2 protein ab54415

製品の詳細

製品名	Recombinant Human FAM19A2 protein
精製度	> 95 % SDS-PAGE. Purity : Greater than 98% by SDS-PAGE gel and HPLC analyses. Endotoxin Level : < 0.1 ng per µg (1EU/µg).
エンドトキシン・レベル	< 0.100 Eu/µg
発現系	Escherichia coli
アクセッション番号	<u>Q8N3H0</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	ANHHKAHHVKTGTCEVVALHRCCNKNKIEERSQTVKCSCFPG QVAGTTTRA APSCVDASIVEQKWWCHMQPCLEGECKVL PDRKGWSCSSGN KVKTTRVT H
予測される分子量	11 kDa
領域	31 to 131

特性

Our **Abpromise guarantee** covers the use of **ab54415** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Lyophilized
備考	Previously labelled as TAFE 2.

前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
再構成	For lot specific reconstitution information please contact our Scientific Support Team.

関連情報

組織特異性	Brain-specific.
配列類似性	Belongs to the FAM19/TAFA family.
細胞内局在	Cytoplasm.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors