

Recombinant human FKBP2 protein ab93681

2 References [画像数 1](#)

製品の詳細

製品名	Recombinant human FKBP2 protein
生理活性	Specific activity is > 700nmol/min/mg, and is defined as the amount of enzyme that cleaves 1nmol of suc-AAPF-pNA per minute at 37°C in Tris-HCl pH 8.0 using chymotrypsin.
精製度	> 90 % SDS-PAGE. ab93681 is purified using conventional chromatography techniques.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MATGAEGKRKLQIGVKRVDHCPIKSRKGDVLHMHYTGKLED GTEFDSSL PQNQPFVFSLGTGQVIKGWDQGLLGMCEGEKRKLVIPSELGY GERGAPPK IPGGATLVFEVELLKIERRTEL

特性

Our **Abpromise guarantee** covers the use of **ab93681** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
	Functional Studies

製品の状態	Liquid
-------	--------


前処理および保存

保存方法および安定性	Shipped on Dry Ice. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: 0.0154% DTT, 0.242% Tris, 10% Glycerol (glycerin, glycerine) This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------	---

関連情報

機能	PPases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides.
組織特異性	T-cells and thymus.
配列類似性	Belongs to the FKBP-type PPase family. FKBP2 subfamily. Contains 1 PPase FKBP-type domain.
細胞内局在	Endoplasmic reticulum membrane.

画像



15% SDS-PAGE analysis of ab93681 (3μg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors