

Recombinant Human GAL4 protein ab89234

画像数 1

製品の詳細

製品名	Recombinant Human GAL4 protein
精製度	> 95 % SDS-PAGE. ab89234 is purified using conventional chromatography techniques.
エンドキシン・レベル	< 1.000 Eu/μg
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human

特性

Our **Abpromise guarantee** covers the use of **ab89234** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot
	SDS-PAGE

製品の状態	Liquid
-------	--------

前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. pH: 8.00 Constituents: 0.0154% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine)
------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------


関連情報

機能	Galectin that binds lactose and a related range of sugars. May be involved in the assembly of adherens junctions.
配列類似性	Contains 2 galectin domains.

ドメイン

Contains two homologous but distinct carbohydrate-binding domains.

画像



15% SDS-PAGE showing ab89234 at approximately 38.1kDa (3μg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors