

# Recombinant Human Galectin 13/PLAC8 protein ab124309

画像数 1

### 製品の詳細

製品名	Recombinant Human Galectin 13/PLAC8 protein
精製度	> 90 % SDS-PAGE. ab124309 was purified using Ni-NTA chromatography. The product is filtered (0.4 µm) but not sterile, hence needs to be filtered by an appropriate sterile filter before using in cell culture.
発現系	Escherichia coli
アクセッション番号	<u><b>Q9UHV8</b></u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	<div> <div>MKHHHHHHAS</div> <div>IIKGTPIHSF</div> <div>AFRFRVHFGN</div> <div>VPFEDGKQFE</div> <div>GFVHRIPPSF</div> </div> <div> <div>MSSLPVPYKL</div> <div>INDPQLQVDF</div> <div>HVVMNRREFG</div> <div>LCIYVHYNEY</div> <div>VKMVQVSRDI</div> </div> <div> <div>PVSLSVGSCV</div> <div>YTDMEEDSDI</div> <div>IWMLEETTDY</div> <div>EIKVNGIRIY</div> <div>SLTSVCVCN</div> </div>

### 特性

Our **Abpromise guarantee** covers the use of **ab124309** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot
	ELISA
	SDS-PAGE
製品の状態	Lyophilized
備考	This product was previously labelled as Galectin 13

## 前処理および保存

### 保存方法および安定性

Shipped at 4°C. Store at -80°C.

Constituents: 99% Phosphate Buffer, 0.44% Sodium chloride

### 再構成

Add deionized water to prepare a working stock solution of 0.5 mg/mL and let the lyophilized pellet dissolve completely. Avoid storage of lyophilized product at room temperature. Aliquot reconstituted protein to avoid repeated freezing/thawing cycles and store at -80°C for long term. Reconstituted protein can be stored at 4°C for a week.

## 関連情報

### 機能

Has lysophospholipase activity.


### 組織特異性

Detected in adult and fetal spleen, fetal kidney, adult urinary bladder and placenta.

### 配列類似性

Contains 1 galectin domain.

## 画像



14% SDS-PAGE separation of Human Galectin 13/PLAC8 protein (ab124309)

Lane 1: reduced and boiled sample (5 µg)

Lane 2: non-reduced and non-boiled sample (5 µg)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors