

### Recombinant Human GC1q R protein ab78834

画像数 1

#### 製品の詳細

製品名	Recombinant Human GC1q R protein
精製度	> 95 % SDS-PAGE. ab78834 is purified using conventional chromatography techniques.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MLHTDGDKAF VDFLSDEIKE ERKIQKHKTL PKMSGGWELE LNGTEAKLVR KVAGEKITVT FNINNSIPT FDGEEEPSQG QKVEEQEPEL TSTPNFVVEV IKNDDGKKAL VLDCHYPEDE VGQEDAEASD IFSIREVSFQ STGESEWKDT NYTLNTDSL D WLYDHLMD F LADRGVDNTF ADELVELSTA LEHQEYITFL EDLKS FVKSQ

#### 特性

Our **Abpromise guarantee** covers the use of **ab78834** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション SDS-PAGE

製品の状態 Liquid


#### 前処理および保存

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.  
pH: 7.50  
Constituents: 0.0154% DTT, 0.242% Tris, 20% Glycerol (glycerin, glycerine)

#### 関連情報

機能	Binds to the globular "heads" of C1Q thus inhibiting C1 activation.
配列類似性	Belongs to the MAM33 family.
細胞内局在	Mitochondrion matrix. Nucleus. Might also be nuclear in some cell types.

## 画像



15% SDS-PAGE showing ab78834 (3μg).

Note: Real molecular weight on SDS-PAGE is shifted up compared to predicted molecular weight.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors