

Recombinant Human GGCX protein ab114658

画像数 1

製品の詳細

製品名	Recombinant Human GGCX protein
発現系	Wheat germ
アクセッション番号	P38435
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	ADFPGLHLENFVSEDLGNTSIQLLQGEVTVELVAEQKNQTLR EGEKMQLP AGEYHKVYTTSPSPSCYMYVYVNTTELALAEQDLAYLQELKEK VENGs
予測される分子量	36 kDa including tags
領域	533 to 629

特性

Our **Abpromise guarantee** covers the use of **ab114658** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot
	SDS-PAGE
	ELISA
製品の状態	Liquid


前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
------------	---

関連情報

機能	Mediates the vitamin K-dependent carboxylation of glutamate residues to calcium-binding gamma-carboxyglutamate (Gla) residues with the concomitant conversion of the reduced hydroquinone form of vitamin K to vitamin K epoxide.
関連疾患	Defects in GGCX are a cause of combined deficiency of vitamin K-dependent clotting factors type 1 (VKCFD1) [MIM:277450]; also known as multiple coagulation factor deficiency III (MCFD3). VKCFD leads to a bleeding tendency that is usually reversed by oral administration of vitamin K. Defects in GGCX are the cause of pseudoxanthoma elasticum-like disorder with multiple coagulation factor deficiency (PXEL-MCFD) [MIM:610842]. This syndrome is characterized by hyperlaxity of the skin involving the entire body. Important phenotypic differences with classical PXE include much more severe skin laxity with spreading toward the trunk and limbs with thick, leathery skin folds rather than confinement to flexural areas, and no decrease in visual acuity. Moreover, detailed electron microscopic analyzes revealed that alterations of elastic fibers as well as their mineralization are slightly different from those in classic PXE.
配列類似性	Belongs to the vitamin K-dependent gamma-carboxylase family.
細胞内局在	Endoplasmic reticulum membrane.

画像



SDS-PAGE analysis of ab114658 on a 12.5% gel stained with Coomassie Blue.

SDS-PAGE - Recombinant Human GGCX protein
(ab114658)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors