

# Recombinant Human GMFB protein ab54143

## 1 References

### 製品の詳細

製品名	Recombinant Human GMFB protein
精製度	> 95 % SDS-PAGE. Purity : Greater than 98% by SDS-PAGE gel and HPLC analyses.
エンドトキシン・レベル	< 0.100 Eu/μg
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	SESLVVCDDVA EDLVEKLRKF RFRKETNNAA IIMKIDKDKR LVLDEELEG ISPDELKDEL PERQPRFIVY SYKYQHDDGR VSYPLCFIFS SPVGCKPEQQ MMYAGSKNKL VQTAELTKVF EIRNTEDLTE EWLREKLGFF H

### 特性

Our **Abpromise guarantee** covers the use of **ab54143** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Lyophilized
備考	This product was previously labelled as GMF beta

### 前処理および保存

保存方法および安定性	Upon delivery aliquot. Upon reconstitution add a carrier protein (0.1% BSA). Store at -20°C. Aliquot the reconstituted solution and re-lyophilize. The re-lyophilized aliquoted peptide should be stored dry at 4°C, preferably in a desiccator. For long term storage it is recommended to add a carrier protein on reconstitution (0.1% HSA or BSA). Working aliquots stored with a carrier protein are stable for at least 3 months at -20°C to -80°C..
------------	--

## 再構成

For lot specific reconstitution information please contact our Scientific Support Team.

## 関連情報

---

機能	This protein causes differentiation of brain cells, stimulation of neural regeneration, and inhibition of proliferation of tumor cells.
配列類似性	Belongs to the actin-binding proteins ADF family. GMF subfamily. Contains 1 ADF-H domain.
翻訳後修飾	Phosphorylated; stimulated by phorbol ester.

---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors