

Recombinant Human Grp75/MOT protein ab115711

画像数 2

製品の詳細

製品名	Recombinant Human Grp75/MOT protein
精製度	>= 80 % SDS-PAGE. ab115711 is purified by multi-step chromatography
発現系	Escherichia coli
アクセッション番号	P38646
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
予測される分子量	75 kDa
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab115711** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot SDS-PAGE
製品の状態	Liquid
備考	Previously labelled as Grp75.

前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. Constituents: 0.79% Tris HCl, 3% Glycerol (glycerin, glycerine), 0.88% Sodium chloride
------------	---

関連情報

機能	Implicated in the control of cell proliferation and cellular aging. May also act as a chaperone.
----	--

配列類似性


Belongs to the heat shock protein 70 family.

翻訳後修飾

Phosphorylated upon DNA damage, probably by ATM or ATR.

細胞内局在

Mitochondrion.

画像

SDS-PAGE analysis:


Lane 1: MW marker

Lane 2: ab115711 at 1 μ g

Lane 3: ab115711 at 2 μ g

Lane 4: ab115711 at 5 μ g

SDS-PAGE - Recombinant Human Grp75/MOT protein (ab115711)



All lanes : Grp75 / MOT antibody

Lane 1 : Recombinant Human Grp75/MOT protein (ab115711) at 0.1 μ g

Lane 2 : HeLa cell lysate (heat shocked)

Western blot - Recombinant Human Grp75/MOT protein (ab115711)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors