

# Recombinant human GST3 / GST pi protein ab86703

画像数 2

### 製品の詳細

製品名	Recombinant human GST3 / GST pi protein
生理活性	Specific activity is > 80 units/mg, and is defined as the amount of enzyme that conjugate 1.0 umole of 1-chloro-2,4-dinitrobenzene (CDNB) with reduced glutathione per minute at pH 6.5 at 25C.
精製度	> 95 % SDS-PAGE. ab86703 is purified using conventional chromatography techniques. Endotoxin Level: < 1.0 EU per 1 µg of protein (determined by LAL method).
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MRGSHHHHHH GMASMTGGGQ MGRDLYDDDD KDRWGSHMPP YTVVYFPVRG RCAALRMLLA DQQQSWKEEV VTVETWQEGS LKASCLYGQL PKFQDGD LTL YQSNTILRHL GRTLGLYGKD QQEAALVDMV NDGVEDLRCK YISLIYTNYE AGKDDYVKAL PGQLKPFETL LSQNQGKTF IVGDQISFAD YNLLDLLLLIH EVLAPGCLDA FPLLSAYVGR LSARPKLKAF LASPEYVNLP INGNKGKQ

### 特性

Our **Abpromise guarantee** covers the use of **ab86703** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Western blot
製品の状態	Liquid

### 前処理および保存

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw

cycles.

pH: 7.00


Constituents: 0.316% Tris HCl, 30% Glycerol (glycerin, glycerine), 0.03% EDTA, 0.002% PMSF

This product is an active protein and may elicit a biological response in vivo, handle with caution.

## 関連情報

機能	Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.
配列類似性	Belongs to the GST superfamily. Pi family. Contains 1 GST C-terminal domain. Contains 1 GST N-terminal domain.

## 画像



15% SDS-PAGE showing ab86703 at approximately 27.6kDa (3μg).

SDS-PAGE - Recombinant Human GST3 / GST pi protein (ab86703)



Western blot - Recombinant Human GST3 / GST pi protein (ab86703)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors