

Recombinant Human HDAC3 protein ab82071

1 References [画像数 1](#)

製品の詳細

製品名	Recombinant Human HDAC3 protein
生理活性	1 unit equals 1 nanogram of purified protein.
精製度	> 95 % SDS-PAGE. Affinity purified
発現系	Baculovirus infected insect cells
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human

特性

Our **Abpromise guarantee** covers the use of **ab82071** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid
備考	1 unit equals 1 nanogram of purified protein.

前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.9 Constituents: 0.75% Potassium chloride, 0.0154% DTT, 0.316% Tris HCl, 0.00584% EDTA, 20% Glycerol (glycerin, glycerine)
------------	---

関連情報

機能	Responsible for the deacetylation of lysine residues on the N-terminal part of the core histones (H2A, H2B, H3 and H4). Histone deacetylation gives a tag for epigenetic repression and plays an important role in transcriptional regulation, cell cycle progression and developmental events.
----	---

Histone deacetylases act via the formation of large multiprotein complexes. Probably participates in the regulation of transcription through its binding to the zinc-finger transcription factor YY1; increases YY1 repression activity. Required to repress transcription of the POU1F1 transcription factor. Acts as a molecular chaperone for shuttling phosphorylated NR2C1 to PML bodies for sumoylation.

組織特異性

Widely expressed.

配列類似性

Belongs to the histone deacetylase family, HD type 1 subfamily.


翻訳後修飾

Sumoylated in vitro.

細胞内局在

Nucleus.

画像



Standard Curve for HDAC3 (Analyte: **HDAC3 protein (Tagged) (ab82071)**); dilution range 1 pg/ml to 1 ug/ml using Capture Antibody **Mouse monoclonal [7G6C5] to HDAC3 (ab76359)** at 5 ug/ml and Detector Antibody **Rabbit polyclonal to HDAC3 (ab63353)** at 0.5 ug/ml.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors