

# Recombinant Human Heme Oxygenase 1 protein ab86919

### 画像数 1

#### 製品の詳細

製品名	Recombinant Human Heme Oxygenase 1 protein		
精製度	> 95 % SDS-PAGE.		
発現系	Escherichia coli		
タンパク質長	Protein fragment		
Animal free	No		
由来	Recombinant		
生物種	Human		
配列	MERPQPHSMP	QDLSEALKEA	TKEVHTQAEN
	AEFMRNFQKG	QVTRDGFKLV	MASLYHIYVA
	LEEEIERNKE	SPVFAPVYFP	EELHRKAALE
	QDLAFWYGPR	WQEVIPYTPA	MQRYVKRLHE
	VGRTEPELLV	AHAYTRYLGD	LSGGQVLKKI
	AQKALDLPSS	GEGLAFFTFP	NIASATKFKQ
	LYRSRMNSLE	MTPAVRQRVI	EEAKTAFLLN
	IQLFEELQEL	LTHDTKDQSP	SRAPGLRQRA
	SNKVQDSAPV	ETPRGKPLLN	TRSQAPLEHH HHHH
領域	1 to 266		

#### 特性

Our **Abpromise guarantee** covers the use of **ab86919** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid

#### 前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
	pH: 8.00
	Constituents: 0.00174% PMSF, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.29%

## 関連情報


---

機能	Heme oxygenase cleaves the heme ring at the alpha methene bridge to form biliverdin. Biliverdin is subsequently converted to bilirubin by biliverdin reductase. Under physiological conditions, the activity of heme oxygenase is highest in the spleen, where senescent erythrocytes are sequestered and destroyed.
配列類似性	Belongs to the heme oxygenase family.
細胞内局在	Microsome. Endoplasmic reticulum.

---

## 画像

---



showing ab86919 at 31.4 kDa on 15% SDS-PAGE (4µg)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors