

Recombinant Human HINT1 protein ab87362

★★★★★ [1 Abreviews](#) [1 References](#) [画像数 1](#)

製品の詳細

製品名	Recombinant Human HINT1 protein
精製度	> 90 % SDS-PAGE. Purified by using conventional chromatography.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MADEIAKAQV ARPGGDTIFG KIIRKEIPAK IIFEDDRCLA FHDISPQAPT HFLVIPKKHI SQISVAEDDD ESLLGHLMIV GKKCAADLGL NKGYRMVVNE GSDGGQSVYH VHLHVLGGRQ MHWPPG

特性

Our **[Abpromise guarantee](#)** covers the use of **ab87362** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション SDS-PAGE

製品の状態 Liquid

前処理および保存


保存方法および安定性 Shipped at 4°C. Store at +4°C.
pH: 8.00
Constituents: 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine)

関連情報

機能 Hydrolyzes adenosine 5'-monophosphoramidate substrates such as AMP-morpholidate, AMP-N-alanine methyl ester, AMP-alpha-acetyl lysine methyl ester and AMP-NH2.

組織特異性	Widely expressed.
配列類似性	Belongs to the HINT family. Contains 1 HIT domain.
ドメイン	The histidine triad, also called HIT motif, forms part of the binding loop for the alpha-phosphate of purine mononucleotide.
細胞内局在	Cytoplasm. Nucleus. Interaction with CDK7 leads to a more nuclear localization.

画像



ab87362 on 15% SDS-PAGE (3μg)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors