

Recombinant Human HSBP1 protein ab51233

画像数 1

製品の詳細

製品名	Recombinant Human HSBP1 protein
精製度	> 95 % SDS-PAGE. This protein was purified to apparent homogeneity by using conventional column chromatography techniques.
発現系	Escherichia coli
アクセッション番号	<u>NP_001528.1</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MAETDPKTVQ DLTSVVQTLL QQMQDKFQTM SDQIIGRIDD MSSRIDDLK NIADLMTQAG VEELESENKI PATQKS

特性

Our **Abpromise guarantee** covers the use of **ab51233** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid

前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.0154% DTT, 0.316% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.29% Sodium chloride
------------	--

関連情報

機能	Negative regulator of the heat shock response. Negatively affects HSF1 DNA-binding activity.
----	--

May have a role in the suppression of the activation of the stress response during the aging process.


配列類似性

Belongs to the HSBP1 family.

細胞内局在

Nucleus.

画像



ab51233 run on a 15% SDS-APGE gel with molecular weight markers.

SDS-PAGE - Recombinant Human HSBP1 protein
(ab51233)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors