

Recombinant Human HSPBAP1 protein ab107135

画像数 1

製品の詳細

製品名	Recombinant Human HSPBAP1 protein
精製度	> 80 % SDS-PAGE. ab107135 is purified using conventional chromatography techniques.
発現系	Escherichia coli
アクセッション番号	<u>Q96EW2</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHSSGLVPRGSHMAAGSEATTPVIVAAGAGGEEG EHVKPFKP EKAKEIIMSLQQPAIFCNMVFDPARHWNAYLSQVLHGKQI RFRMGMK MSTVPQFETTCNYVEATLEEFLLTWNCDSISGPFDRDYDHSK FWAYADYK YFVSLFEDKTDLFQDVKWSDFGFPGRNGQESTLWIGSLGAHT PCHLDSYG CNLVFQVQGRKRWHLPPEDETPFLYPTRIPIYEESVFSKINV VNPDLKRF PQFRKAQRHAVTLSPGQVLFVPRHWHYVESIDPVTVSINSW IELEEDHL ARVEEAITRMLVCALKTAENPQNTRAWLNPTEVEETSHAVNC CYLNAAVS AFFDRCRTSEVVEIQALRTDGEHMKKEELNVCNHMEVGTGS QNLTTGTD KPEAASPFGLPVAQRSEEPSEGGIFGSDGKDFVDDKG EHFGKLHC AKRQQIMSNSENAIEEQIASNTTTTPQTFISTDDLLDCLVNP QVTRIVAQLLIQGRSL
予測される分子量	57 kDa
領域	1 to 488
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab107135** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
	Mass Spectrometry
質量分析	MALDI-TOF
製品の状態	Liquid


前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
	pH: 7.40
	Constituent: PBS

関連情報

機能	May play a role in cellular stress response.
組織特異性	Widely expressed.
関連疾患	Note=A chromosomal aberration involving HSPBAP1 has been found in a family with renal carcinoma. Translocation t(2;3)(q35;q21) with the putative pseudogene DIRC3. Produces an hybrid mRNA encoding a truncated HSPBAP1 lacking the first 36 amino acids.
配列類似性	Contains 1 JmjC domain.
細胞内局在	Cytoplasm.

画像



15% SDS-PAGE showing ab107135 (3µg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors