

Recombinant human IL-20 protein ab49819

製品の詳細

製品名	Recombinant human IL-20 protein
生理活性	Biological activity determined by its ability to activate STAT following receptor ligand interaction.
精製度	> 98 % SDS-PAGE.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MLKTLNLGSC VIATNLQEIR NGFSEIRGSV QAKDGNIDIR ILRRTESLQD TKPANRCCLL RHLLRLYLDR VFKNYQTPDH YTLRKISSLA NSFLTIKKDL RLCHAHMTCH CGEEAMKKYS QILSHFEKLE PQAAVVKALG ELDILLQWME ETE

特性

Our **Abpromise guarantee** covers the use of **ab49819** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Functional Studies
製品の状態	Lyophilized

前処理および保存

保存方法および安定性	Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C long term. Constituent: 0.16% Sodium phosphate This product is an active protein and may elicit a biological response in vivo, handle with caution.
再構成	Reconstitute to 1mg/ml using sterile water.

関連情報

機能	Cytokine that may be involved in epidermal function and psoriasis. Acts through STAT3.
組織特異性	Expressed at very low levels in skin, trachea, and other tissues.
配列類似性	Belongs to the IL-10 family.
細胞内局在	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors