

### Recombinant Human Inhibin beta A protein ab119449

#### 製品の詳細

製品名	Recombinant Human Inhibin beta A protein
精製度	> 95 % SDS-PAGE. Purity is >95% by SDS PAGE and HPLC analysis.
発現系	Escherichia coli
アクセッション番号	<b>P08476</b>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
予測される分子量	14 kDa including tags
領域	311 to 426
タグ	His tag N-Terminus

#### 特性

Our **Abpromise guarantee** covers the use of **ab119449** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot ELISA SDS-PAGE HPLC
製品の状態	Liquid
備考	ab119449 should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

#### 前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Constituents: 49% Tris buffered saline, 50% Glycerol
------------	---

## 関連情報

---

機能	Inhibins and activins inhibit and activate, respectively, the secretion of follitropin by the pituitary gland. Inhibins/activins are involved in regulating a number of diverse functions such as hypothalamic and pituitary hormone secretion, gonadal hormone secretion, germ cell development and maturation, erythroid differentiation, insulin secretion, nerve cell survival, embryonic axial development or bone growth, depending on their subunit composition. Inhibins appear to oppose the functions of activins.
配列類似性	Belongs to the TGF-beta family.
細胞内局在	Secreted.

---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors