

### Recombinant Human Inhibin beta A protein ab53506

画像数 1

#### 製品の詳細

製品名	Recombinant Human Inhibin beta A protein
精製度	> 90 % SDS-PAGE.
発現系	Escherichia coli
アクセッション番号	<b>P08476</b>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
領域	311 to 426
タグ	His tag N-Terminus
配列の追加情報	The sequence is the mature full length without the signal peptide and propeptide. Hexahistidine tag at N-terminus.

#### 特性

Our **Abpromise guarantee** covers the use of **ab53506** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid

#### 前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
------------	---

#### 関連情報

機能	Inhibins and activins inhibit and activate, respectively, the secretion of follitropin by the pituitary gland. Inhibins/activins are involved in regulating a number of diverse functions such as hypothalamic and pituitary hormone secretion, gonadal hormone secretion, germ cell development and maturation, erythroid differentiation, insulin secretion, nerve cell survival,
----	---

embryonic axial development or bone growth, depending on their subunit composition. Inhibins appear to oppose the functions of activins.


**配列類似性**

Belongs to the TGF-beta family.

**細胞内局在**

Secreted.

**画像**



**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors