

Recombinant Human IRAKM protein ab125546

画像数 1

製品の詳細

製品名	Recombinant Human IRAKM protein
精製度	> 80 % SDS-PAGE. Assessed by densitometry. Affinity purified.
発現系	Baculovirus infected Sf9 cells
アクセッション番号	Q9Y616
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
予測される分子量	96 kDa
領域	1 to 596

特性

Our **Abpromise guarantee** covers the use of **ab125546** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot SDS-PAGE
----------	--------------------------

製品の状態	Liquid
-------	--------

前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.31% Glutathione, 0.002% PMSF, 0.004% DTT, 0.79% Tris HCl, 0.003% EDTA, 25% Glycerol (glycerin, glycerine), 0.29% Sodium chloride
------------	---

関連情報

機能	Inhibits dissociation of IRAK1 and IRAK4 from the Toll-like receptor signaling complex by either
----	--

inhibiting the phosphorylation of IRAK1 and IRAK4 or stabilizing the receptor complex.

組織特異性

Expressed predominantly in peripheral blood lymphocytes.

関連疾患

Defects in IRAK3 are associated with susceptibility to asthma-related traits type 5 (ASRT5) [MIM:611064]. Asthma-related traits include clinical symptoms of asthma, such as coughing, wheezing, dyspnea, bronchial hyperresponsiveness as assessed by methacholine challenge test, serum IgE levels, atopy, and atopic dermatitis.


配列類似性

Belongs to the protein kinase superfamily. TKL Ser/Thr protein kinase family. Pelle subfamily.

Contains 1 death domain.

Contains 1 protein kinase domain.

画像



SDS-PAGE analysis of ab125546.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors