

Recombinant human IRR protein ab69988

画像数 2

製品の詳細

製品名	Recombinant human IRR protein
生理活性	Specific activity: 111 nmol/min/mg as per kinase assay.
精製度	> 90 % SDS-PAGE. Affinity purified.
発現系	Baculovirus infected Sf9 cells
アクセッション番号	<u>P14616</u>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
予測される分子量	65 kDa
領域	945 to 1297
タグ	GST tag N-Terminus
配列の追加情報	Recombinant human fragment of Insulin Receptor R (945-end).

特性

Our **Abpromise guarantee** covers the use of **ab69988** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Functional Studies SDS-PAGE
製品の状態	Liquid
備考	<u>ab204853</u> (IRS1 peptide) can be utilized as a substrate for assessing kinase activity

前処理および保存


保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.0038% EGTA, 0.00174% PMSF, 0.00385% DTT, 0.79% Tris HCl, 0.00292% EDTA, 25% Glycerol (glycerin, glycerine), 0.87% Sodium chloride
------------	--

This product is an active protein and may elicit a biological response in vivo, handle with caution.

関連情報


機能	This receptor probably binds an insulin related protein and has a tyrosine-protein kinase activity. It phosphorylates the insulin receptor substrates IRS-1 and IRS-2.
配列類似性	Belongs to the protein kinase superfamily. Tyr protein kinase family. Insulin receptor subfamily. Contains 3 fibronectin type-III domains. Contains 1 protein kinase domain.
細胞内局在	Membrane.

画像



IRR enzymatic activity.

Functional Studies - Recombinant human IRR protein (ab69988)



SDS page analysis of ab69988

SDS-PAGE - Recombinant human IRR protein (ab69988)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors