

Product datasheet

Recombinant Human LETMD1/HCCR-1 protein ab79566

画像数 1

製品の詳細

製品名	Recombinant Human LETMD1/HCCR-1 protein
精製度	> 90 % SDS-PAGE.
発現系	Escherichia coli
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
領域	183 to 330

特性

Our **Abpromise guarantee** covers the use of **ab79566** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot ELISA SDS-PAGE
製品の状態	Liquid
備考	This product was previously labelled as LETMD1

前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Constituent: 2.68% PBS
------------	-----------------------------------------------------------------------------------------------------------------

関連情報


関連性	LETMD1 is involved in tumorigenesis and may function as a negative regulator of the TP53/p53. It
-----	--------------------------------------------------------------------------------------------------

is found in kidney, liver, skeletal muscle, heart and brain and is over-expressed in various tumors including leukemia, lymphoma, and carcinomas of the breast, kidney, ovary, stomach, colon and uterine cervix. There are six isoforms.

細胞内局在

Membrane; Single-pass membrane protein

画像



ab79566 on 15% SDS PAGE gel, Coomassie blue staining.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors