

Recombinant Human LHPP protein ab116175

[2 References](#) [画像数 1](#)

製品の詳細

製品名	Recombinant Human LHPP protein		
精製度	> 90 % SDS-PAGE. ab116175 was purified using conventional chromatography.		
発現系	Escherichia coli		
アクセッション番号	Q9H008		
タンパク質長	Full length protein		
Animal free	No		
由来	Recombinant		
生物種	Human		
配列	MRGSHHHHHP	WYASMTGGQQ	MGRDLYDDDD
	KDRWGSHMAP	WGKRLAGVRG	VLLDISGVLY
	DSGAGGGTAI	AGSVEAVARL	KRSRLKVRFC
	TNESQKSRAE	LVGQLQRLGF	DISEQEVTP
	APAACQILKE	QGLRPYLLIH	DGVRSEFDQI
	DTSNPNCVVI	ADAGESFSYQ	NMNAFQVLM
	ELEKPVLISL	GKGYYKETS	GLMLDVGPM
	KALEYACGIK	AEVVGKPSPE	FFKSALQAIG
	VEAHQAVMIG	DDIVGDVGG	QRCGMRALQV
	RTGKFRPSDE	HHPEVKADGY	VDNLAEAVDL LLQHADK
予測される分子量	34 kDa including tags		
領域	1 to 270		
タグ	His tag N-Terminus , DDDDK tag N-Terminus		

特性

Our **Abpromise guarantee** covers the use of **ab116175** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
	Mass Spectrometry
質量分析	MALDI-TOF
製品の状態	Liquid

前処理および保存

保存方法および安定性

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.02% (R*,R*)-1,4-Dimercaptobutan-2,3-diol, 0.32% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride, 0.002% PMSF

関連情報


組織特異性

Expressed in liver, kidney and moderately in brain.

配列類似性

Belongs to the HAD-like hydrolase superfamily.

画像



ab116175 on a 15% SDS-PAGE (3ug)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors