

### Recombinant Human LI Cadherin protein ab114173

画像数 1

#### 製品の詳細

製品名	Recombinant Human LI Cadherin protein
発現系	Wheat germ
アクセッション番号	<b>Q12864</b>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	EGKFSGPLKPMTFSIYEGQEPSQIIFQFKANPPAVTFELTGE TDNIFVIE REGLLYYNRALDRETRSTHNLQVAALDANGIIVEGPVPITIE VKDINDNR PTFLQSKY
予測される分子量	38 kDa including tags
領域	24 to 131

#### 特性

Our **Abpromise guarantee** covers the use of **ab114173** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
	Western blot
	ELISA
製品の状態	Liquid


#### 前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
------------	---

#### 関連情報

<b>機能</b>	Cadherins are calcium dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types. LI-cadherin may have a role in the morphological organization of liver and intestine. Involved in intestinal peptide transport.
<b>組織特異性</b>	Expressed in the gastrointestinal tract and pancreatic duct. Not detected in kidney, lung, liver, brain, adrenal gland and skin.
<b>配列類似性</b>	Contains 7 cadherin domains.
<b>細胞内局在</b>	Cell membrane.

## 画像



ab114173 on a 12.5% SDS-PAGE Stained with Coomassie Blue.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors