

Recombinant Human LSM1 protein ab101157

画像数 1

製品の詳細

| | |
|-------------|--|
| 製品名 | Recombinant Human LSM1 protein |
| 精製度 | > 95 % SDS-PAGE. ab101157 is purified using conventional chromatography techniques. |
| 発現系 | Escherichia coli |
| アクセッション番号 | <u>O15116</u> |
| タンパク質長 | Full length protein |
| Animal free | No |
| 由来 | Recombinant |
| 生物種 | Human |
| 配列 | MGSSHHHHHHSSGLVPRGSHMNYPGTASLIEDIDKKHLVLL RDGRTLIG FLRSIDQFANLVLHQTVERIHVGKKYGDI PRGIFVVRGENVV LLGEIDLE KESDTPLQQVSIEEILEEQRVEQQTKEAEKLVQALKDRGL SIPRADTL DEY |
| 予測される分子量 | 17 kDa including tags |
| 領域 | 1 to 133 |
| タグ | His tag N-Terminus |

特性

Our **Abpromise guarantee** covers the use of **ab101157** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|----------|-------------------|
| アプリケーション | SDS-PAGE |
| | Mass Spectrometry |
| 質量分析 | MALDI-TOF |
| 製品の状態 | Liquid |

前処理および保存

| | |
|------------|---|
| 保存方法および安定性 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw |
|------------|---|

cycles.


pH: 7.50

Constituents: 0.0154% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

関連情報

| | |
|-------|--|
| 機能 | Plays a role in replication-dependent histone mRNA degradation. Binds specifically to the 3'-terminal U-tract of U6 snRNA. |
| 組織特異性 | Has elevated expression in pancreatic cancer and in several cancer-derived cell lines. |
| 配列類似性 | Belongs to the snRNP Sm proteins family. |
| 細胞内局在 | Nucleus. |

画像



15% SDS-PAGE showing ab101157 at approximately 17.3kDa (3μg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors