

### Recombinant Human Lymphotactin/ATAC protein ab85406

画像数 1

#### 製品の詳細

製品名	Recombinant Human Lymphotactin/ATAC protein
精製度	> 95 % SDS-PAGE. ab85406 is purified using conventional chromatography techniques.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHH SSGLVPRGSH MVGSEVSDKR TCVSLTTQRL PVSRIKTYTI TEGSLRAVIF ITKRGLKVCA DPQATWVRDV VRSMRKSNT RNNMIQTKPT GTQQSTNTAV TLTG

#### 特性

Our **Abpromise guarantee** covers the use of **ab85406** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
質量分析	MALDI-TOF
製品の状態	Liquid
備考	This product was previously labelled as Lymphotactin


#### 前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.  pH: 8.00 Constituents: 0.0308% DTT, 0.316% Tris HCl, 30% Glycerol (glycerin, glycerine), 1.16% Sodium chloride
------------	--

## 関連情報

機能	Chemotactic activity for lymphocytes but not for monocytes or neutrophils.
組織特異性	Highest level in spleen, lower in peripheral leukocytes and very low levels in lung, colon and small intestine.
配列類似性	Belongs to the intercrine gamma family.
細胞内局在	Secreted.

## 画像



15% SDS-PAGE showing ab85406 (3μg).

Note: Real molecular weight on SDS-PAGE is shifted up compared to predicted molecular weight.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors