

Recombinant Human MDH1 protein ab99244

[1 References](#) [画像数 1](#)

製品の詳細

製品名	Recombinant Human MDH1 protein
精製度	> 95 % SDS-PAGE. ab99244 is purified using conventional chromatography techniques.
エンドキシン・レベル	< 1.000 Eu/ μ g
発現系	Escherichia coli
アクセッション番号	P40925
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MSEPIRVLVTGAAGQIAYSLLYSIGNGSVFGKDQPIILVLLD ITPMMGVL DGVLMELQDCALPLLKDVIATDKEDVAFKDLDVAILVGSMPR REGMERKD LLKANVKIFKSQGAALDKYAKKSVKIVVGNPANTNCLTASK SAPSIPKE NFSCLTRLDHNRKAQIALKLGVTANDVKNVIIWGNHSSTQY PDVNHAKV KLGQKEVGVYEALKDDSWLKGEFVTTVQQRGA AVIKARKLSS AMSAAKAI CDHVRDIWFGTPEGEFVSMGVISDGNSYGV PDDL LYSFPVVI KNKTWKFV EGLPINDFSREKMDLTAKELTEEKESAFEFLSSALEHHHHH H
予測される分子量	37 kDa including tags
領域	1 to 334
タグ	His tag C-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab99244** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
	Mass Spectrometry
質量分析	MALDI-TOF
製品の状態	Liquid


前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. pH: 8.00 Constituents: 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine)
------------	--

関連情報

配列類似性	Belongs to the LDH/MDH superfamily. MDH type 2 family.
翻訳後修飾	ISGylated.
細胞内局在	Cytoplasm.

画像



SDS-PAGE analysis of ab99244

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors