

Recombinant Human MED3 protein (denatured) ab180327

1 References [画像数 1](#)

製品の詳細

製品名	Recombinant Human MED3 protein (denatured)
精製度	> 80 % SDS-PAGE.
発現系	Escherichia coli
アクセッション番号	Q6P2C8
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHHSSGLVPRGSHMGSMADVINVSVNLEAFSQAIS AIQALRSS VSRVFDCLKDGMRNKETLEGREKAFIAHFQDNLHSVNRDLNE LERLSNLV GKPSENHPLHNSGLLSLDPVQDKTPLYSQLLQAYKWSNKLQY HAGLASGL LNQQLKRSANQMGVSAKRRPKAQPTTLVLPPQYVDDVISRI DRMFPEMS IHLSRPNGTSAMLLVTLGKVLKVIVVMRSLFIDRTIVKGYNE NVYTEDGK LDIWSKSNYQVFQKVTDHATTALLHYQLQMPDVVVRSFMTW LRSYIKLF QAPCQRCGKFLQDGLPPTWRDFRTL EAFHDTCRQ
予測される分子量	38 kDa including tags
領域	1 to 311
タグ	His tag N-Terminus
配列の追加情報	NP_004260

特性

Our **Abpromise guarantee** covers the use of **ab180327** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
----------	----------

製品の状態	Liquid
備考	This product was previously labelled as CRSP8
前処理および保存	
保存方法および安定性	<p>Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.</p> <p>pH: 8.00</p> <p>Constituents: 10% Glycerol (glycerin, glycerine), 2.4% Urea, 0.32% Tris HCl</p>
関連情報	
機能	<p>Component of the Mediator complex, a coactivator involved in the regulated transcription of nearly all RNA polymerase II-dependent genes. Mediator functions as a bridge to convey information from gene-specific regulatory proteins to the basal RNA polymerase II transcription machinery. Mediator is recruited to promoters by direct interactions with regulatory proteins and serves as a scaffold for the assembly of a functional preinitiation complex with RNA polymerase II and the general transcription factors.</p>
配列類似性	Belongs to the Mediator complex subunit 27 family.
細胞内局在	Nucleus.

画像



15% SDS-PAGE analysis of 3µg ab180327.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors