

# Recombinant Human NME4 protein ab107133

画像数 1

### 製品の詳細

---

製品名	Recombinant Human NME4 protein
精製度	> 90 % SDS-PAGE. ab107133 is purified using conventional chromatography techniques.
発現系	Escherichia coli
アクセッション番号	<b><u>O00746</u></b>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHSSGLVPRGSHMPSWTRERTLVAVKPDGVQRRL VGDVIQRF ERRGFTLVGMKMLQAPESVLAEHYQDLRRKPFYPALIRYMSS GPVVAMVW EGYNVVRASRAMIGHTDSAEAAPGTIRGDFSVHISRNVIHAS DSVEGAQR EIQWLFQSSSELVSWADGGQHSSIHPA
予測される分子量	20 kDa including tags
領域	33 to 187
タグ	His tag N-Terminus

### 特性

---

Our **Abpromise guarantee** covers the use of **ab107133** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Mass Spectrometry
質量分析	MALDI-TOF
製品の状態	Liquid

### 前処理および保存

---

保存方法および安定性 Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -

80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.316% Tris HCl, 40% Glycerol (glycerin, glycerine), 1.16% Sodium chloride

## 関連情報


---

機能	Major role in the synthesis of nucleoside triphosphates other than ATP.
組織特異性	Widely distributed. Found at very high levels in prostate, heart, liver, small intestine, and skeletal muscle tissues, and in low amounts in the brain and in blood leukocytes.
配列類似性	Belongs to the NDK family.
細胞内局在	Mitochondrion intermembrane space.

---

## 画像

---



15% SDS-PAGE showing ab107133 (3µg).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors