

Recombinant Human NNT1 protein (denatured) ab111634

画像数 1

製品の詳細

| | |
|-------------|--|
| 製品名 | Recombinant Human NNT1 protein (denatured) |
| 精製度 | > 90 % SDS-PAGE. |
| 発現系 | Escherichia coli |
| アクセッション番号 | <u>Q9UBD9</u> |
| タンパク質長 | Full length protein |
| Animal free | No |
| 由来 | Recombinant |
| 生物種 | Human |
| 配列 | MGSSHHHHHHSSGLVPRGSHMLNRTGDPGPGPSIQKTYDLTR YLEHQLRS LAGTYLNYLGPPFNEPDFNPPRLGAETLPRATVDLEVWRSLN DKLRLTQN YEAYSHLLCYLRGLNRQAATAELRRSLAHFCTSLQGLLGSIA GVMAALGY PLPQPLPGTEPTWTPGPAHSDFLQKMDDFWLLKELQTLWRS AKDFNRLK KKMQPAAAVTLHLGAHGF |
| 予測される分子量 | 25 kDa including tags |
| 領域 | 28 to 225 |
| タグ | His tag N-Terminus |
| 製品の詳細 | Recombinant Human NNT1 protein |

特性

Our **Abpromise guarantee** covers the use of **ab111634** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|----------|----------|
| アプリケーション | SDS-PAGE |
| 製品の状態 | Liquid |

前処理および保存

保存方法および安定性

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 3.50

Constituents: 0.59% Sodium citrate, 2.4% Urea, 10% Glycerol (glycerin, glycerine)

関連情報

機能

Cytokine with B-cell stimulating capability. Binds to and activates the ILST/gp130 receptor.

組織特異性

Expressed predominantly in lymph nodes, spleen, peripheral blood lymphocytes, bone marrow, and fetal liver.

関連疾患

Defects in CLCF1 are the cause of cold-induced sweating syndrome type 2 (CISS2) [MIM:610313]. Cold-induced sweating syndrome (CISS) is an autosomal recessive disorder characterized by profuse sweating induced by cool surroundings (temperatures of 7 to 18 degrees Celsius). Additional abnormalities include a high-arched palate, nasal voice, depressed nasal bridge, inability to fully extend the elbows and kyphoscoliosis.


配列類似性

Belongs to the IL-6 superfamily.

細胞内局在

Secreted.

画像



15% SDS-PAGE showing ab111634 at approximately 24.6kDa (3μg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors