

Recombinant Human NUCB2 protein ab101042

画像数 1

製品の詳細

製品名	Recombinant Human NUCB2 protein
精製度	> 95 % SDS-PAGE. Purity as determined by densitometric image analysis: > 95%
発現系	Escherichia coli
アクセッション番号	P80303
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	MKHHHHHHASVPIDIDKTKVQNIHPVESAKIEPPDTGLYYD EYLKQVIDV LETDKHFREKLQKADIEEIKSGRLSKELDLVSHHVRTKLDEL
予測される分子量	11 kDa including tags
領域	25 to 106
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab101042** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot SDS-PAGE
製品の状態	Lyophilized
備考	Store lyophilized protein at -20°C. Lyophilized protein remains stable until the expiry date when stored at -20°C. Previously labelled as Nucleobindin 2.

前処理および保存

保存方法および安定性 Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -80°C.

Avoid freeze / thaw cycle. Please see notes section.

Constituents: 0.242% Tris, 0.29% Sodium chloride

再構成

Add deionized water to prepare a working stock solution of approximately 0.5 mg/mL and let the lyophilized pellet dissolve completely. Product is not sterile! Please filter the product by an appropriate sterile filter before using it in the cell culture.

関連情報

機能 Calcium-binding protein. May have a role in calcium homeostasis.


組織特異性 Predominantly expressed in spleen, testis and normal stomach.

配列類似性 Belongs to the nucleobindin family.

Contains 2 EF-hand domains.

細胞内局在 Membrane. Cytoplasm. Secreted.

画像



14% SDS-PAGE analysis of ab101042

Lane 1: Molecular weight markers; 14, 21, 31, 45, 66, 97 kDa

Lane 2: ab101042 (reduced and heated, 5µg/lane)

Lane 3: ab101042 (non-reduced and non-heated, 5µg/lane)

SDS-PAGE - Recombinant Human NUCB2 protein
(ab101042)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors