

## Product datasheet

# Recombinant Human p23 protein ab113183

## 1 References [画像数 2](#)

### 製品の詳細

---

製品名	Recombinant Human p23 protein
精製度	> 90 % SDS-PAGE. Purified by multi step chromatography
発現系	Escherichia coli
アクセッション番号	<b><u>Q15185</u></b>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
予測される分子量	23 kDa

### 特性

---

Our **Abpromise guarantee** covers the use of **ab113183** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Western blot
----------	--------------------------

製品の状態	Liquid
-------	--------

### 前処理および保存

---

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. Constituent: 100% PBS
------------	--

### 関連情報

---

機能	Molecular chaperone that localizes to genomic response elements in a hormone-dependent manner and disrupts receptor-mediated transcriptional activation, by promoting disassembly of transcriptional regulatory complexes.
パスウェイ	Lipid metabolism; prostaglandin biosynthesis.

**配列類似性**


Belongs to the p23/wos2 family.

Contains 1 CS domain.

**細胞内局在**

Cytoplasm.


**画像**



Western blot - Recombinant Human p23 protein (ab113183)

Recombinant Human p23 protein (ab113183) at 0.1  $\mu$ g + Antibody which recognises p23

Developed using the ECL technique.



SDS-PAGE - Recombinant Human p23 protein (ab113183)

ab113183, at 0.5 $\mu$ g in lane 1, 1 $\mu$ g in lane 2, 2 $\mu$ g in lane 3 and 5 $\mu$ g in lane 4, staining Human p23 using SDS-PAGE. Imperial stain used.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors