

Recombinant Human P4HA1 protein ab114624

画像数 1

製品の詳細

製品名	Recombinant Human P4HA1 protein
発現系	Wheat germ
アクセッション番号	P13674
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MIWYILIIIGILLPQSLAHPGFFTSIGQMTDLIHTEKDLVTSL KDYIKAEE DKLEQIKKWAEKLDRLTSTATKDPEGFVGHPVNAFKLMKRLN TEWSELEN LVLKDMSDGFISNLTIQRQYFPNDEDQVGAALKLLRLQDTYN LDTDTISK GNLPGVKHKSFLTAEDCFELGKVAYTEADYYHTELWMEQALR QLDEGEIS TIDKVSVDYLSYAVYQQGDLKALLLTKKLELDPEHQAN GNLKYFEY IMAKEKDVNKSASDDQSDQKTPKKKGVAVDYLPERQKYEML CRGEGIKM TPRRQKKLFCRYHDGNRNPKFILAPAKQEDEWDKPRIIRFHD IISDAEIE IVKDLAKPRLRRATISNPITGDLETVHYRISKSAWLSGYENP VVSRIINMR IQDLTGLDVSTAEELQVANYGVGGQYEPHFDFARKDEPDAFK ELGTGNRI ATWLFYMSDVSAGGATVFPEVGASVWPKKGTAVFWYNLFASG EGDYSTRH AACPVLVGNKWSNKWLHERGQEFRRPCTLSELE
予測される分子量	85 kDa including tags
領域	1 to 534

特性

Our **Abpromise guarantee** covers the use of **ab114624** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
	Western blot
	ELISA
製品の状態	Liquid


前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
------------	---

関連情報

機能	Catalyzes the post-translational formation of 4-hydroxyproline in -Xaa-Pro-Gly- sequences in collagens and other proteins.
配列類似性	Belongs to the P4HA family. Contains 1 Fe2OG dioxygenase domain. Contains 1 TPR repeat.
細胞内局在	Endoplasmic reticulum lumen.

画像



ab114624 analysed on a 12.5% SDS-PAGE stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors