

Recombinant Human Pancreatic Polypeptide protein ab114376

画像数 1

製品の詳細

製品名	Recombinant Human Pancreatic Polypeptide protein
発現系	Wheat germ
アクセッション番号	P01298
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MAAARLCLSLLLLSTCVALLLQPLLGAQGAPLEPVYPGDNAT PEQMAQYA ADLRRYINMLTRPRYGKRHKEDTLAFSEWGSPHAAVPRELSP LDL
予測される分子量	36 kDa including tags
領域	1 to 95

特性

Our **Abpromise guarantee** covers the use of **ab114376** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA
	SDS-PAGE
	Western blot
製品の状態	Liquid

備考


前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
------------	---

関連情報

関連性	Pancreatic hormone is synthesized in pancreatic islets of Langerhans and acts as a regulator of pancreatic and gastrointestinal functions and may be important in the regulation of food intake. Plasma levels of this hormone have been shown to be reduced in conditions associated with increased food intake and elevated in anorexia nervosa. In addition, infusion of this hormone in obese rodents has shown to decrease weight gain.
細胞内局在	Secreted

画像



12.5% SDS-PAGE with ab114376 at 36.19kDa, stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors