

## Product datasheet

# Recombinant Human PAR1/Thrombin Receptor protein ab114346

画像数 1

### 製品の詳細

製品名	Recombinant Human PAR1/Thrombin Receptor protein
発現系	Wheat germ
アクセッション番号	<b>P25116</b>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	SFLLRNPNDKYEPFWEDEEKNESGLTEYRLVSINKSSPLQKQ LPAFISED ASGYLTSSWLT
予測される分子量	32 kDa including tags
領域	42 to 102

### 特性

Our **Abpromise guarantee** covers the use of **ab114346** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
	Western blot
	ELISA

製品の状態	Liquid
-------	--------


### 前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
------------	---

### 関連情報

<b>機能</b>	High affinity receptor for activated thrombin coupled to G proteins that stimulate phosphoinositide hydrolysis. May play a role in platelets activation and in vascular development.
<b>組織特異性</b>	Platelets and vascular endothelial cells.
<b>配列類似性</b>	Belongs to the G-protein coupled receptor 1 family.
<b>翻訳後修飾</b>	A proteolytic cleavage generates a new N-terminus that functions as a tethered ligand. Phosphorylated; probably mediating desensitization prior to the uncoupling and internalization of the receptor.
<b>細胞内局在</b>	Cell membrane.

## 画像



ab114346 analysed on a 12.5% SDS-PAGE gel stained with Coomassie Blue.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors