

Recombinant human PF4 protein ab9562

1 References

製品の詳細

製品名	Recombinant human PF4 protein
精製度	> 98 % SDS-PAGE. Endotoxin level is less than 0.1 ng per g (1EU/g).
発現系	Escherichia coli
アクセッション番号	P02776
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	EAEEDGDLQC LCVKTTTSQVR PRHITSLEVI KAGPHCPTAQ LIATLKNGRK ICLDLQAPLY KKIIKKLLES
予測される分子量	8 kDa
領域	31 to 101

特性

Our **Abpromise guarantee** covers the use of **ab9562** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Functional Studies SDS-PAGE
製品の状態	Lyophilized
備考	The biological activity of this product is determined by its ability to chemoattract human fibroblasts using a concentration range of 1.0-10.0 ng/ml.

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------	--

再構成

Reconstitute in 200ul sterile water to a conc. of 0.1 mg/ml.

関連情報

機能	Released during platelet aggregation. Neutralizes the anticoagulant effect of heparin because it binds more strongly to heparin than to the chondroitin-4-sulfate chains of the carrier molecule. Chemotactic for neutrophils and monocytes. Inhibits endothelial cell proliferation, the short form is a more potent inhibitor than the longer form.
配列類似性	Belongs to the intercrine alpha (chemokine CxC) family.
翻訳後修飾	Binds non-covalently to a proteoglycan molecule.
細胞内局在	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors