

Recombinant Human phospholipase A2 protein ab63843

画像数 1

製品の詳細

製品名	Recombinant Human phospholipase A2 protein
精製度	> 95 % SDS-PAGE.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	<p>MRGSHHHHHH GMASHMNLVN FHRMIKLTG</p> <p>KEAALSYGFY GCHCGVGGRG SPKDATDRCC</p> <p>VTHDCCYKRL EKRGCGTKFL SYKFSNSGSR</p> <p>ITCAKQDSCR SQLCECDKAA ATCFARNKTT</p> <p>YNKKYQYYSN KHCRCSTPRC</p>
領域	21 to 144

特性

Our **Abpromise guarantee** covers the use of **ab63843** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot
製品の状態	Lyophilized
備考	Product is NOT sterile! Please filter the product by an appropriate sterile filter before using it in cell culture.

前処理および保存


保存方法および安定性	Shipped at 4°C. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. Constituent: 0.41% Sodium acetate
再構成	Add 0.1M Acetate buffer pH4 to prepare a working stock solution of approximately 0.5 mg/ml and let the lyophilized pellet dissolve completely. For conversion into higher pH value, we recommend intensive dilution by relevant buffer to a concentration of 10µg/ml. In higher concentrations the

solubility of this antigen is limited.

関連情報

機能	Thought to participate in the regulation of the phospholipid metabolism in biomembranes including eicosanoid biosynthesis. Catalyzes the calcium-dependent hydrolysis of the 2-acyl groups in 3-sn-phosphoglycerides.
配列類似性	Belongs to the phospholipase A2 family.
細胞内局在	Membrane.

画像



SDS-PAGE - Recombinant Human phospholipase A2 protein (ab63843)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors