

Recombinant Human Phospholipase D1 protein ab114568

画像数 1

製品の詳細

製品名	Recombinant Human Phospholipase D1 protein
発現系	Wheat germ
アクセッション番号	<u>Q13393</u>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	GYLDDPSEDIQDPVSDKFFKEVWVSTAARNATIYDKVFRCLP NDEVHNL I QLRDFINKPVLAKEDPIRAEEELKKIRGFLVQFPFYFLSEES LLPSVGTK EAIVPMEVWT
予測される分子量	38 kDa including tags
領域	965 to 1074

特性

Our **Abpromise guarantee** covers the use of **ab114568** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
	Western blot
	ELISA
製品の状態	Liquid


前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
------------	---

関連情報

機能	Implicated as a critical step in numerous cellular pathways, including signal transduction, membrane trafficking, and the regulation of mitosis. May be involved in the regulation of perinuclear intravesicular membrane traffic.
組織特異性	Expressed abundantly in the pancreas and heart and at high levels in brain, placenta, spleen, uterus and small intestine.
配列類似性	Belongs to the phospholipase D family. Contains 1 PH domain. Contains 2 PLD phosphodiesterase domains. Contains 1 PX (phox homology) domain.
細胞内局在	Cytoplasm > perinuclear region. Endoplasmic reticulum membrane. Golgi apparatus membrane. Late endosome membrane.

画像



12.5% SDS-PAGE Stained with Coomassie Blue

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors