

# Recombinant Human PI 3 Kinase p55 gamma protein ab127108

画像数 1

### 製品の詳細

---

製品名	Recombinant Human PI 3 Kinase p55 gamma protein
精製度	> 80 % SDS-PAGE. Purity was determined by densitometry. Affinity purified.
発現系	Baculovirus infected Sf9 cells
アクセッション番号	<b><u>Q8N381</u></b>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
予測される分子量	61 kDa including tags
領域	1 to 461
タグ	His tag N-Terminus

### 特性

---

Our **Abpromise guarantee** covers the use of **ab127108** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot SDS-PAGE
製品の状態	Liquid

### 前処理および保存

---


保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.00 Preservative: 1.02% Imidazole Constituents: 0.002% PMSF, 0.07% Sodium phosphate, 0.004% DTT, 25% Glycerol (glycerin, glycerine), 1.76% Sodium chloride
------------	--

### 関連情報

---

<b>機能</b>	Binds to activated (phosphorylated) protein-tyrosine kinases through its SH2 domain and regulates their kinase activity. During insulin stimulation, it also binds to IRS-1.
<b>組織特異性</b>	Highest levels in brain and testis. Lower levels in adipose tissue, kidney, heart, lung and skeletal muscle.
<b>配列類似性</b>	Belongs to the PI3K p85 subunit family. Contains 2 SH2 domains.

## 画像



SDS-PAGE analysis of ab127108.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors