

## Product datasheet

# Recombinant Human Placental alkaline phosphatase (PLAP) protein ab114268

画像数 1

### 製品の詳細

製品名	Recombinant Human Placental alkaline phosphatase (PLAP) protein
発現系	Wheat germ
アクセッション番号	<b><u>P05187</u></b>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MLGPCMLLLLLLLGLRLQLSLGIIPVEEENPDFWNREAAEAL GAAKKLQP AQTAAKNLIIIFLDGMGVSTVTAARILKGQKKDKLGPEIPLA MDRFPYVA LSKTYNVDKHVPDSGATATAYLCGVKGNFQTIGLSAAARFNQ CNTTRGNE VISVMNRAKKAGKSVGVTTRVQHASPAGTYAHTVNRNWYS DADVPASA RQEGCQDIATQLISNMDIDVILGGGRKYMFRMGTPDPEYPDD YSQGGTRL DGKNLVQEWLAKRQGARYVWNRTLMQASLDPSVTHLMGLFE PGDMKYEI HRDSTLDPSLMEMTEAALRLLSRNPRGFFLFVEGGRIDHGH ESRAYRAL TETIMFDDAIERAGQLTSEEDTSLVLTADHSHVFSFGGYPLR GSSIFGLA PGKARDRKAYTVLLYNGPGYVLKDGARPDVTESESGSPEYR QQSAVPLD EETHAGEDVAVFARGPQAHLVHGVQEQTFFIAHVMAFAACLEP YTACDLAP PAGTTDAAHPGRSVVPALLPLLAGTLLLLLETATAP
予測される分子量	85 kDa including tags
領域	1 to 535
タグ	GST tag N-Terminus

## 特性

---

Our **Abpromise guarantee** covers the use of **ab114268** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション                    ELISA  
    SDS-PAGE  
    Western blot

製品の状態                        Liquid

## 前処理および保存

---

保存方法および安定性            Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

## 関連情報

---


配列類似性                        Belongs to the alkaline phosphatase family.

細胞内局在                        Cell membrane.

---

## 画像

---



12.5% SDS-PAGE analysis of Human Placental alkaline phosphatase (PLAP) protein (ab114268). Stained with Coomassie Blue.

SDS-PAGE - Recombinant Human Placental alkaline phosphatase (PLAP) protein (ab114268)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors