

# Recombinant Human PRDM14 protein ab196396

[1 References](#) [画像数 1](#)

### 製品の詳細

製品名	Recombinant Human PRDM14 protein
精製度	>= 70 % SDS-PAGE.
発現系	Escherichia coli
アクセッション番号	<b><u>Q9GZV8</u></b>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	MSPILGYWKIKGLVQPTRLLEYLEEKYEEHLYERDEGDKWR NKKFELGL EFPNLPYYIDGDVKLTQSMAIIRYIADKHMLGGCPKERA EI SMLEGAVL DIRYGVSRIAYS KDFETLKVDFLSKLPEMLKMFEDRLCHKTY LNGDHVTH PDFMLYDALDVVLYMDPMCLDAFPKLVCFKKRIEAI PQIDKY LKSSKYIA WPLQGWQATFGGGDHPKSDLEVL FQ^GPLGSALPRPSEAVP QDKVCYPP ESSPQNLAAYYTPFPSYGHYRNSLATVEEDFQPFRQLEAAAS AAPAMPPF PFRMAPPLLS PGLGLQREPLYDLPWYSKLPWYPIPHVPREV PPFLSSSH EYAGASSED LGHQIIGGDNESGPCCGPD TLIPPPADASLLP EGLRTSQL LPCSPSKQSEDGPKPSNQEKGSPARFQFTEEDLHFVLYGVTP SLEHPASL HHAISGLLVPPDSSGSDSLPQTLDKDSLQLPEGLCLMQTVFG EVPHFGVF CSSFIAKGVRF GPFQKVVNASEVKTYGDNSVMWEIFEDGHL SHFIDGKG GTGNWMSYVNCARFPKEQNLVAVQCQGHIFYESCKEIHQNQE LLVWYGDC YEKFLDIPVSLQVTEPGKQPSGPSEESAEGYRCERC GKVFTY

予測される分子量	72 kDa including tags
領域	2 to 412
タグ	GST tag N-Terminus

### 特性

Our **Abpromise guarantee** covers the use of **ab196396** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid


### 前処理および保存

保存方法および安定性	Shipped on Dry Ice. Store at -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: 0.63% Tris HCl, 0.64% Sodium chloride, 0.02% Potassium chloride, 0.04% Tween, 20% Glycerol (glycerin, glycerine)
------------	--

### 関連情報

機能	May be involved in transcriptional regulation.
配列類似性	Contains 6 C2H2-type zinc fingers. Contains 1 SET domain.
細胞内局在	Nucleus.

### 画像



4-20% SDS-PAGE with Coomassie staining: 2 µg ab196396.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors