

Recombinant human PTP alpha protein ab61138

画像数 5

製品の詳細

製品名	Recombinant human PTP alpha protein
精製度	> 90 % Densitometry.
発現系	Escherichia coli
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
領域	174 to 802

特性

Our **Abpromise guarantee** covers the use of **ab61138** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
	Functional Studies

製品の状態	Liquid
-------	--------


前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.0038% EGTA, 0.00174% PMSF, 0.00385% DTT, 0.79% Tris HCl, 0.00292% EDTA, 25% Glycerol (glycerin, glycerine), 0.87% Sodium chloride This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------	--

関連情報


配列類似性	Belongs to the protein-tyrosine phosphatase family. Receptor class 4 subfamily. Contains 2 tyrosine-protein phosphatase domains.
-------	---

画像



Functional Studies - Recombinant human PTP alpha protein (ab61138)


The specific activity of PTP alpha (ab61138) was determined to be 2117 nmol/min/mg as per activity assay protocol



SDS-PAGE - Recombinant human PTP alpha protein (ab61138)


SDS PAGE analysis of ab61138

SDS PAGE analysis of ab61138




SDS-PAGE - Recombinant human PTP alpha protein (ab61138)

Sample Kinase Activity Plot.



Functional Studies - Recombinant human PTP alpha protein (ab61138)

ab61138 on SDS-PAGE, MW ~96 kDa.



SDS-PAGE - Recombinant human PTP alpha protein (ab61138)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors