

Recombinant human PTPRS protein ab61149

画像数 2

製品の詳細

製品名	Recombinant human PTPRS protein
精製度	> 90 % Densitometry.
発現系	Escherichia coli
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
予測される分子量	63 kDa including tags
領域	883 to 1210
タグ	GST tag N-Terminus
配列の追加情報	GeneBank Accession number: NM_130855

特性

Our **Abpromise guarantee** covers the use of **ab61149** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Functional Studies
製品の状態	Liquid


前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.50 Constituents: 0.418% MOPS, 0.00174% PMSF, 0.00385% DTT, 30% Glycerol (glycerin, glycerine), 0.29% Sodium chloride This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------	--

関連情報

機能	Interacts with LAR-interacting protein LIP.1.
組織特異性	Detected in all tissues tested except for placenta and liver.
配列類似性	Belongs to the protein-tyrosine phosphatase family. Receptor class 2A subfamily. Contains 8 fibronectin type-III domains. Contains 3 Ig-like C2-type (immunoglobulin-like) domains. Contains 2 tyrosine-protein phosphatase domains.
翻訳後修飾	A cleavage occurs, separating the extracellular domain from the transmembrane segment. This process called 'ectodomain shedding' is thought to be involved in receptor desensitization, signal transduction and/or membrane localization.
細胞内局在	Membrane.

画像



Sample Kinase Activity Plot.



ab61149 on SDS-PAGE, MW ~63 kDa.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors