

Recombinant human Rab2 protein ab90797

画像数 1

製品の詳細

製品名	Recombinant human Rab2 protein
生理活性	Activity: GTPgammaS binding: >750 mmol/mol.
精製度	> 90 % SDS-PAGE.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human

特性

Our **Abpromise guarantee** covers the use of **ab90797** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Functional Studies
	SDS-PAGE
	Western blot

製品の状態	Liquid
-------	--------

前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.20 Constituents: 0.443% GDP, 0.019% Magnesium chloride, 0.077% DTT, 0.595% HEPES, 0.232% Sodium chloride This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------	--

関連情報

機能	Required for protein transport from the endoplasmic reticulum to the Golgi complex.
----	---


配列類似性

Belongs to the small GTPase superfamily. Rab family.

細胞内局在

Endoplasmic reticulum-Golgi intermediate compartment membrane. Melanosome. Endoplasmic reticulum membrane. Golgi apparatus membrane. Identified by mass spectrometry in melanosome fractions from stage I to stage IV.

画像



Anti-Rab2 antibody ([ab96509](#)) at 1/500 dilution + Recombinant human Rab2 protein (ab90797) at 0.01 µg

Secondary

Goat Anti-Rabbit IgG H&L (HRP) preadsorbed ([ab97080](#)) at 1/5000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

Exposure time: 20 seconds

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors