

Recombinant Human RABIF protein ab100886

画像数 1

製品の詳細

製品名	Recombinant Human RABIF protein
精製度	> 85 % SDS-PAGE. ab100886 is purified using conventional chromatography techniques.
発現系	Escherichia coli
アクセッション番号	<u>P47224</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHH SSGLVPRGSH MEPAEQPSELVSAEGRNRKA VLCQRCGSRV LQPGTALFSR RQLFLPSMRK KPALSDG SNP DGDLLQEHWL VEDMFIFENV GFTKDVGN IK FLVCADCEIG PIGWHCLDDK NSFYVALERV SHE
予測される分子量	16 kDa including tags
領域	1 to 123
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab100886** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Mass Spectrometry
	SDS-PAGE
質量分析	MALDI-TOF
製品の状態	Liquid

前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
------------	---

pH: 8.00


Constituents: 0.0154% DTT, 0.316% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

関連情報

関連性

RABIF belongs to the DSS4/MSS4 family. It is a guanine-nucleotide-releasing protein that acts on members of the SEC4/YPT1/RAB subfamily. It stimulates GDP release from both YPT1 and RAB3A, but is less active on these proteins than on the SEC4 protein. RABIF might play a general role in vesicular transport. Ubiquitously expressed.

画像



ab100886 (3 µg) on an SDS-PAGE gel (15 %).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

-
- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors