

### Recombinant Human Ran (mutated T24N) protein ab90206

#### 製品の詳細

製品名	Recombinant Human Ran (mutated T24N) protein
精製度	> 95 % SDS-PAGE. Protein preparation is 18% GDP-loaded, measured by HPLC.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
修飾	mutated T24N

#### 特性

Our **Abpromise guarantee** covers the use of **ab90206** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid
備考	Ran T24N mutation results in no or little binding to GTP and is known to act as a negative inhibitor of nuclear import.

#### 前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.00 Constituents: 0.1001% DTE (1,4-Dithioerythritol), 1.125% Potassium chloride, 0.42% Potassium phosphate, 0.0475% Magnesium chloride, 10% Glycerol
------------	---

#### 関連情報

機能	GTP-binding protein involved in nucleocytoplasmic transport. Required for the import of protein into the nucleus and also for RNA export. Involved in chromatin condensation and control of cell
----	--

	cycle (By similarity). The complex with BIRC5/ survivin plays a role in mitotic spindle formation by serving as a physical scaffold to help deliver the RAN effector molecule TPX2 to microtubules. Enhances AR-mediated transactivation. Transactivation decreases as the poly-Gln length within AR increases.
<b>組織特異性</b>	Expressed in a variety of tissues.
<b>配列類似性</b>	Belongs to the small GTPase superfamily. Ran family.
<b>翻訳後修飾</b>	The N-terminus is blocked.
<b>細胞内局在</b>	Nucleus. Cytoplasm. Melanosome. Becomes dispersed throughout the cytoplasm during mitosis. Identified by mass spectrometry in melanosome fractions from stage I to stage IV.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors